

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Klebsiella pneumoniae subsp. pneumoniae Autoinducer 2 import ATP-binding protein LsrA(IsrA),Escherichia Coli Preis auf Anfrage CSB-EP412572KAX

|                            |                                                                                                                                                                                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Recombinant Klebsiella pneumoniae subsp. pneumoniae Autoinducer 2 import ATP-binding protein LsrA(IsrA),Escherichia Coli Preis auf Anfrage                                                                                                                                    |
| Biozol Catalog Number      | CSB-EP412572KAX                                                                                                                                                                                                                                                               |
| Supplier Catalog Number    | CSB-EP412572KAX                                                                                                                                                                                                                                                               |
| Alternative Catalog Number | CSB-EP412572KAX-1,CSB-EP412572KAX-100,CSB-EP412572KAX-20                                                                                                                                                                                                                      |
| Manufacturer               | Cusabio                                                                                                                                                                                                                                                                       |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                                              |
| Species Reactivity         | Bacteria                                                                                                                                                                                                                                                                      |
| Tag                        | The tag type will be determined during production process.                                                                                                                                                                                                                    |
| UniProt                    | <a href="#">A6TEB8</a>                                                                                                                                                                                                                                                        |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                   |
| Source                     | E.coli                                                                                                                                                                                                                                                                        |
| Purity                     | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                   |
| Form                       | Lyophilized powder                                                                                                                                                                                                                                                            |
| Sequence                   | MMKPLLEARQ ICKQFSGVAV LKGIDFTLGA GQVHALMGGN GAGKSTLMKI<br>IAGVETPD SG ELTIGERAFA RLNPALAHQL GIYLV PQEPM LFPNLSVREN<br>ILFRLPKRAD TAARLQDKLQ QLNCQINLDA SASTLEVADQ QMVEILRGLM<br>REARILILDE PTASLTPGET ERLFSQIRAL QALDVGIVFI SHKLPEIRQL<br>ASHISVMRDG AVVLSGETAA YRDEQLINAM TP |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|