

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant *Xenopus tropicalis* Coiled-coil domain-containing protein 160 homolog, *Escherichia Coli* Preis auf Anfrage CSB-EP429231XBF

|                            |                                                                                                                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Recombinant <i>Xenopus tropicalis</i> Coiled-coil domain-containing protein 160 homolog, <i>Escherichia Coli</i> Preis auf Anfrage                                                                                                                                           |
| Biozol Catalog Number      | CSB-EP429231XBF                                                                                                                                                                                                                                                              |
| Supplier Catalog Number    | CSB-EP429231XBF                                                                                                                                                                                                                                                              |
| Alternative Catalog Number | CSB-EP429231XBF-1, CSB-EP429231XBF-100, CSB-EP429231XBF-20                                                                                                                                                                                                                   |
| Manufacturer               | Cusabio                                                                                                                                                                                                                                                                      |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                                             |
| Species Reactivity         | Frog                                                                                                                                                                                                                                                                         |
| Tag                        | The tag type will be determined during production process.                                                                                                                                                                                                                   |
| UniProt                    | <a href="#">A9JSR5</a>                                                                                                                                                                                                                                                       |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                  |
| Source                     | E.coli                                                                                                                                                                                                                                                                       |
| Purity                     | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                  |
| Form                       | Lyophilized powder                                                                                                                                                                                                                                                           |
| Sequence                   | MEDNKHWIEK LFPPHFTVQD FFSQSFEPEQ LISEKISKEK VKNVEKMYHN<br>AAERCLEAGK WKRKEELSNI FSYDPNIPDD QKEKEPPVGN LHDNVTCPA<br>AHVDDQCIWR ENELKVL RHE MQQRYSEGA K FKNQLDACKL<br>ELSELKAKDK NTEQELGKAK EALTLSKTHI RNKGILVKQL QKDKLQKDSE<br>IQSLKKDLHE KSVMINSLNK NLCIAGEEIQ ELKLKSKDLE QE |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|