

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Polynucleobacter sp. Phosphatidylserine decarboxylase proenzyme(psd),Escherichia Coli Preis auf Anfrage CSB-EP018033PYU

| | |
|----------------------------|--|
| Article Name | Recombinant Polynucleobacter sp. Phosphatidylserine decarboxylase proenzyme(psd),Escherichia Coli Preis auf Anfrage |
| Biozol Catalog Number | CSB-EP018033PYU |
| Supplier Catalog Number | CSB-EP018033PYU |
| Alternative Catalog Number | CSB-EP018033PYU-1, CSB-EP018033PYU-100, CSB-EP018033PYU-20 |
| Manufacturer | Cusabio |
| Category | Proteine/Peptide |
| Tag | The tag type will be determined during production process. |
| UniProt | A4SXR2 |
| Buffer | Tris/PBS-based buffer, 6% Trehalose, pH 8.0 |
| Source | E.coli |
| Purity | Greater than 85% as determined by SDS-PAGE. |
| Form | Lyophilized powder |
| Sequence | MMYPHPPIAK EGWPYLALIG IVTLPVHYIG GIAWSWPLWI IFIFVLQFFR DPQRIPAMGR DLVLSPADGR IVVVEKTNDP YADREALKIS VFMNVFNVHS NRSVNGLVK EIQYFPGKFV NADLDKASTE NERNAVVIDA NGHIVTLVQV AGLIARRILC YIHVGDRLLKA GERYGFIRFG |

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.