

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Danio rerio RING finger protein 141(rnf141),Escherichia Coli Preis auf Anfrage CSB-EP019839DIL

Article Name	Recombinant Danio rerio RING finger protein 141(rnf141),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP019839DIL
Supplier Catalog Number	CSB-EP019839DIL
Alternative Catalog Number	CSB-EP019839DIL-1, CSB-EP019839DIL-100, CSB-EP019839DIL-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Fish
Tag	The tag type will be determined during production process.
UniProt	Q6IV56
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MGQQLSGQAV TRLPEKLIKH VGLVRDSGYL TYDEFLGRVA ELNDVTAKLA SGQKKHLLFE VQPGSDSSAL WKVAVRIVCT KINKEDGMVE ASRIMNLYQF IQLYKDITSQ AAEVFSSNVA AEGSSEDTCQ ASMWMGRVKQ LTDEEECCIC MDGKADLILP CAHSFCQKCI DKWSGQSRNC PVCRIQVTAA NESWVMSDAP TGEDVAGYIL NLADEAGHPH RP

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--