

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Human Serpin B10 (SERPINB10) CSB-EP021066HU

Article Name	Recombinant Human Serpin B10 (SERPINB10)
Biozol Catalog Number	CSB-EP021066HU
Supplier Catalog Number	CSB-EP021066HU
Alternative Catalog Number	CSB-EP021066HU-1,CSB-EP021066HU-100,CSB-EP021066HU-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Human
Molecular Weight	52.3 kDa
Tag	C-terminal 6xHis-tagged
UniProt	P48595
Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Source	E.coli
Expression System	1-397aa
Purity	Greater than 90% as determined by SDS-PAGE.
Form	Liquid or Lyophilized powder

Sequence	MDSLATSINQFALELSKKLAESAQGKNIFFSSWSISTSLTIVYLGAKGTTAAQM AQVLQFNRDQGVKCDPESEKKRKMEFNLSNSEEIHSDFQTLISEILKPNDDYLL KTANAIYGEKTYAFHNKYLEDMKTYFGAEPQPVNFVEASDQIRKDINSWVERQ TEGKIQNLLPDDSVDSTTRMILVNALYFKGIWEHQFLVQNTTEKPFRINETTSK PVQMMFMKKKLHIFHIEKPKAVGLQLYYKSRDLSLLILLP
Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.