

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Recombinant Bovine Pulmonary surfactant-associated protein D(SFTPD), Escherichia Coli Preis auf Anfrage CSB-EP021175BO

Article Name	Recombinant Bovine Pulmonary surfactant-associated protein D(SFTPD),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP021175BO
Supplier Catalog Number	CSB-EP021175BO
Alternative Catalog Number	CSB-EP021175BO-1, CSB-EP021175BO-100, CSB-EP021175BO-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bovine
Tag	The tag type will be determined during production process.
UniProt	P35246
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	AEMKIYSQKT MANACTLVMC SPPEDGLPGR DGRDGREGPR GEKGDPGSPG PAGRAGMPGP AGPIGLKGDN GSAGEPGPKG DTGPPGPPGM PGPAGREGPS GKQGSMGPPG TPGPKGDTGP KGGVGAPGIQ GSPGPAGLKG ERGAPGEPGA PGRAGAPGPA GAIGPQGPSG ARGPPGLKGD RGTPGERGAK GESGLAEVNA LRQRVGILEG QLQRLQNAFS QYKKAMLFPN GR

## **Application Notes**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.