

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Bothrops jararaca Platelet-aggregating proteinase PA-BJ, Escherichia Coli Preis auf Anfrage CSB-EP305691BUW

Article Name	Recombinant Bothrops jararaca Platelet-aggregating proteinase PA-BJ, Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP305691BUW
Supplier Catalog Number	CSB-EP305691BUW
Alternative Catalog Number	CSB-EP305691BUW-1, CSB-EP305691BUW-100, CSB-EP305691BUW-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	P81824
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	VVGGR PCKINVHRSL VLLYNSSSL CSGTLINQEW VLTAHCDSK NFKMKLGVHS IKIRNKNERT RHPKEKFICP NRKKDDVLDK DIMLIRLNRP VSNSEHIAPL SLPSSPPSVG SVCYVMGWGK ISSTKETYPD VPHCAKINIL DHAVCRAAYT WWPATSTTLC AGILQGGKDT CEGDSGGPLI CNGLQGIVSG GGNPCGQPRK PALYTKVFDY LPWIESIAG TTTATCP

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.