

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Recombinant Streptococcus pyogenes serotype M3 UPF0473 protein SpyM3\_1796 (SpyM3\_1796),Escherichia Coli Preis auf Anfrage CSB-EP319758SMV

| Article Name               | Recombinant Streptococcus pyogenes serotype M3 UPF0473 protein SpyM3_1796 (SpyM3_1796),Escherichia Coli Preis auf Anfrage |
|----------------------------|---|
| Biozol Catalog Number      | CSB-EP319758SMV   |
| Supplier Catalog Number    | CSB-EP319758SMV   |
| Alternative Catalog Number | CSB-EP319758SMV-1, CSB-EP319758SMV-100, CSB-EP319758SMV-<br>20  |
| Manufacturer               | Cusabio   |
| Category                   | Proteine/Peptide  |
| Tag                        | The tag type will be determined during production process.  |
| UniProt                    | PODH26  |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0   |
| Source                     | E.coli  |
| Purity                     | Greater than 85% as determined by SDS-PAGE.   |
| Form                       | Lyophilized powder  |
| Sequence                   | MTHNHENDHQ HEVITLVDEQ GNETLFEILL TIDGREEFGK NYVLLVPAGS<br>EEDESGEIEI QAYSFTENED GTEGDLQPIP EDSDAEWDMI EEVFNSFLDE N        |

**Application Notes** 

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.