

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

**Recombinant Streptococcus pyogenes serotype M3
Phosphoribosylformylglycinamide cyclo-ligase(purM),Escherichia Coli
Preis auf Anfrage
CSB-EP320467SMV**

Article Name	Recombinant Streptococcus pyogenes serotype M3 Phosphoribosylformylglycinamide cyclo-ligase(purM),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP320467SMV
Supplier Catalog Number	CSB-EP320467SMV
Alternative Catalog Number	CSB-EP320467SMV-1, CSB-EP320467SMV-100, CSB-EP320467SMV-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	P0DD62
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder

Sequence	MSEKNAYAKS GVDVEAGYEV VERIKKHVAR TERAGVMGVL GGFGGMFDLS KTGVPKEPVLV SGTG DGVGTKL MLAIKYDKHD TIGQDCVAMC VNDIIAAGAE PLYFLDYIAT GKNNPVKLEE VVSGVAEGCV QAGVALIGGE TAEMPGMYGE DDYDLAGFAV GVAEKSQIID GSKVKEGDIL LGLASSGIHS NGYSLVRRVF ADYTGKELLP ELEGKQLKDV LLEPTRIYVK AA
Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.