

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant *Vibrio alginolyticus* Trk system potassium uptake protein trkA(trkA),*Escherichia Coli* Preis auf Anfrage CSB-EP334562VEU

Article Name	Recombinant <i>Vibrio alginolyticus</i> Trk system potassium uptake protein trkA(trkA), <i>Escherichia Coli</i> Preis auf Anfrage
Biozol Catalog Number	CSB-EP334562VEU
Supplier Catalog Number	CSB-EP334562VEU
Alternative Catalog Number	CSB-EP334562VEU-1, CSB-EP334562VEU-100, CSB-EP334562VEU-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	P39448
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MKIILGAGQ VGGTLAENLV GENNDITIVD NNADRLRELQ DKYDLRVVNG HASHPDVLHE AGAQDADMLV AVTNTDETNM AACQVAFTLF NTPNRVARIR SPEYLAEKEA LFKSGAIPVD HLIAPEELVT SYIERLIQYP GALQVVSAE QKVS LVAVKA YYGGPLVGNA LSALREHMPH IDTRVAAIFR QGRPIRPQGT TIEADDEVF FVAASNHIRS VMSELQRLEK PY

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--