

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant *Vibrio splendidus* UDP-2,3-diacylglycosamine hydrolase(IpxH),*Escherichia Coli* Preis auf Anfrage CSB-EP492760VFH

Article Name	Recombinant <i>Vibrio splendidus</i> UDP-2,3-diacylglycosamine hydrolase(IpxH), <i>Escherichia Coli</i> Preis auf Anfrage
Biozol Catalog Number	CSB-EP492760VFH
Supplier Catalog Number	CSB-EP492760VFH
Alternative Catalog Number	CSB-EP492760VFH-1, CSB-EP492760VFH-100, CSB-EP492760VFH-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	B7VGZ3
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MKTYFISDLH LAPSRQDITD CFLTFMKNEA LEADALYVLG DLFEFWIGDD DTSEFANSIR QTFIDLVNTG VPCYFTQG NR DFLVGKKFAK QTGVQLLDEV STIDIYGQKA VVLHGDTLCT EDIKYLAFRE KVHQPWLQWV FNRIPFFIKK KIVSKVQSDI KDDKQTKSLD IMDVTQQEVE NVMKQHNVDL MIHGHTHRPN IHSFSANNCT KTRIVLGDWY TQGSVLVFTP QS

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.