

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Desulfovibrio desulfuricans UDP-N-acetylglucosamine 1-carboxyvinyltransferase(murA),Escherichia Coli Preis auf Anfrage CSB-EP493497DJN

Article Name	Recombinant Desulfovibrio desulfuricans UDP-N-acetylglucosamine 1-carboxyvinyltransferase(murA),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP493497DJN
Supplier Catalog Number	CSB-EP493497DJN
Alternative Catalog Number	CSB-EP493497DJN-1, CSB-EP493497DJN-100, CSB-EP493497DJN-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	B8J283
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MDKLVIEGGV PLTGSIEVSG SKNAALPILF AAILPEEPVT ITNVPDLRDI HTTLNLLKVL GCDQCQYENGQ VRIVPGSLLP EAPYDLVRTM RASVLCGLPL LARIGQARVA LPGGCAIGAR PVDQHLKGLE QMGASFQLEE GYIIGRCRKL TGAHITFDMP TVGGTENLLM AAVLAEGKTV LENVALEPEV VDLANFLCAC GARISGQGTS CIRIEGVTSL HQATYPVMPD RI

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.