

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Paracoccus denitrificans Protein moxJ(moxJ),Escherichia Coli Preis auf Anfrage CSB-EP341105PBZ

Article Name	Recombinant Paracoccus denitrificans Protein moxJ(moxJ),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP341105PBZ
Supplier Catalog Number	CSB-EP341105PBZ
Alternative Catalog Number	CSB-EP341105PBZ-1, CSB-EP341105PBZ-100, CSB-EP341105PBZ-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	P29900
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	DTTNLR VCASTKDA PF SDAQGAGFEN KIAQVLADEM GATLDLVMLE KDAIYLVRDG IEKDLCDVLV GVDAGDERLL TTRPYRSGY AFVTRQDRNF EGDKWQDVDQ EGFDTFSYRL HSPAETILKY TGRYEYNLIY QASLTNFEDR RNKYTQVEAS RVITEVADGG ADLAIVFAPE AARYVRDSRE PLRMTLITNE IERSDGVII P LQYSQSVGVS KTHPELLGPI EQALQS

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--