

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Clostridium botulinum DNA integrity scanning protein DisA(disA),Escherichia Coli Preis auf Anfrage CSB-EP501282DUE

Article Name	Recombinant Clostridium botulinum DNA integrity scanning protein DisA(disA),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP501282DUE
Supplier Catalog Number	CSB-EP501282DUE
Alternative Catalog Number	CSB-EP501282DUE-1, CSB-EP501282DUE-100, CSB-EP501282DUE-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Clostridium
Tag	The tag type will be determined during production process.
UniProt	C3KVS9
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MRIEKDKELM NILKIMAPGT PLREGLENIL RAKTGGLLIL GDSDQILKLV DGGFKINSEY SPSYVYELAK MDGSIVLSSD LKKILCANAQ LIPDSSIPTF ETGTRHRTAD RVAKQTGAIV IAISQRRNII TVYKGGIKYV LRDSSIIILAR ANQALQTLEK YVAVLDRVVN NLNILEFKDI ATLFVDLTAI QRSEMVMRIV SEIERYICEL GNEGRLIDMQ LSELIKSVEE DG

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--