

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Nostoc sp. Magnesium-chelatase subunit ChII(chII),Escherichia Coli Preis auf Anfrage CSB-EP348278NHR

Article Name	Recombinant Nostoc sp. Magnesium-chelatase subunit ChII(chII),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP348278NHR
Supplier Catalog Number	CSB-EP348278NHR
Alternative Catalog Number	CSB-EP348278NHR-1, CSB-EP348278NHR-100, CSB-EP348278NHR-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	P58571
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MTPTAQTAS ARRVPFPTA IVGQEEMKLA LLLNVIDPKI GGVMIMGDRG TGKSTTIRAL ADLLPEIPVV ANDPFNSDPS DPLMSDEVR QKSGTGAEIP IEFKKVQMVD LPLGATEDRV CGTIDIEKAL SEGVKAFEPG LLAKANRGIL YVDEVNLLDD HLVDVLLDSA ASGWNTVERE GISIRHPARF VLVGSGNPPEE GELRPQLLDR FGMHAEIHTV KEPALRVQIV EQ

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.