

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Exiguobacterium sp. FMN-dependent NADH-azoreductase(azoR), Escherichia Coli Preis auf Anfrage CSB-EP504653EPU

Article Name	Recombinant Exiguobacterium sp. FMN-dependent NADH-azoreductase(azoR), Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP504653EPU
Supplier Catalog Number	CSB-EP504653EPU
Alternative Catalog Number	CSB-EP504653EPU-1, CSB-EP504653EPU-100, CSB-EP504653EPU-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	C4LOW8
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MTKVLYVSAN PKPTLSYSK QVAETFVSTL KAENASIEVE AIELYDVDVQ EIDGDVLSAW GKFASGEALT DVEAKKVGTM SGMLEKFMEA DLYVFATPMW NFFFPARMKM FLDSVLMAGK TFRYTEQGPGV GLENKQAIH IQGTGGIYTG TDLNFADAYL RQALAFVGVS EVTTVAVEGM NQYPDKIEEI VADAKAKAEA LAKEVAGAVT V

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--