

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Escherichia coli O6 4-hydroxy-3-methylbut-2-enyl diphosphate reductase(ispH),Escherichia Coli Preis auf Anfrage CSB-EP351185EGX

|                            |                                                                                                                                                                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Recombinant Escherichia coli O6 4-hydroxy-3-methylbut-2-enyl diphosphate reductase(ispH),Escherichia Coli Preis auf Anfrage                                                                                                                                                 |
| Biozol Catalog Number      | CSB-EP351185EGX                                                                                                                                                                                                                                                             |
| Supplier Catalog Number    | CSB-EP351185EGX                                                                                                                                                                                                                                                             |
| Alternative Catalog Number | CSB-EP351185EGX-1, CSB-EP351185EGX-100, CSB-EP351185EGX-20                                                                                                                                                                                                                  |
| Manufacturer               | Cusabio                                                                                                                                                                                                                                                                     |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                                            |
| Species Reactivity         | E. coli                                                                                                                                                                                                                                                                     |
| Tag                        | The tag type will be determined during production process.                                                                                                                                                                                                                  |
| UniProt                    | <a href="#">P62624</a>                                                                                                                                                                                                                                                      |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                 |
| Source                     | E.coli                                                                                                                                                                                                                                                                      |
| Purity                     | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                 |
| Form                       | Lyophilized powder                                                                                                                                                                                                                                                          |
| Sequence                   | MQILLANPRG FCAGVDRAIS IVENALAIYG APIYVRHEVV HNRYVVDSLRE<br>ERGAIFIEQI SEVPDGAILI FSAHGVSQAV RNEAKSRDLT VFATCPLVT<br>KVHMEVARAS RRGEESILIG HAGHPEVEGT MGQYSNPEGG<br>MYLVESPDDV WKLTVKNEEK LSFMTQTTLV VDDTSDVIDA LRKRFPKIVG<br>PRKDDICYAT TNRQEAVRAL AEQAEVVLVV GSKNSSNSNR LA |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|