

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Arthroderma gypseum Neutral protease 2 homolog MGYG_03465 (MGYG_03465), Escherichia Coli Preis auf Anfrage CSB-EP516755DOO

Article Name	Recombinant Arthroderma gypseum Neutral protease 2 homolog MGYG_03465 (MGYG_03465), Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP516755DOO
Supplier Catalog Number	CSB-EP516755DOO
Alternative Catalog Number	CSB-EP516755DOO-1, CSB-EP516755DOO-100, CSB-EP516755DOO-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	E4US45
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	TK IVSGSCNQNT LRATMNALGN SARLAQAASR AASQNPRKFQ EYFRTNDANA KQRVIARLNS VARESSSANA GVTYYCSDT MGGCKPRVLA YTLPSRNLVV NCPIYYNLPP LTKQCHAQDQ ATTTLHEFTH NPAVASPHCQ DYAYGYQQCI SLPAAKAVQN ADNYALFANG MLSNLFVFTL N

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.