

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Shigella flexneri Iron uptake system component EfeO(efeO),Escherichia Coli Preis auf Anfrage CSB-EP364596SZB

Article Name	Recombinant Shigella flexneri Iron uptake system component EfeO(efeO),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP364596SZB
Supplier Catalog Number	CSB-EP364596SZB
Alternative Catalog Number	CSB-EP364596SZB-1, CSB-EP364596SZB-100, CSB-EP364596SZB-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	P0AB25
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	ADVP QVKVTVTDKQ CEPMTITVNA GKTQFIIQNH SQKALEWEIL KGV MVVEERE NIAPGFSQKM TANLQPG EYD MTCGLLTNPK GKLIVKGEAT ADAAQSDALL SLGGAITAYK AYVMAETTQL VTDTKAFTDA IKAGDIEKAK ALYAPTRQHY ERIEPIAELF SLDLGSIDAR EDDYEQKAAD PKFTGFHRLE KALFGDNTTK GMDQYAEQLY TDVVDLQKRI SELAFPPS

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--