

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Arthroderma benhamiae Methionine aminopeptidase 2 homolog ARB_02691 (ARB_02691), Escherichia Coli Preis auf Anfrage CSB-EP520564DOM

Article Name	Recombinant Arthroderma benhamiae Methionine aminopeptidase 2 homolog ARB_02691 (ARB_02691), Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP520564DOM
Supplier Catalog Number	CSB-EP520564DOM
Alternative Catalog Number	CSB-EP520564DOM-1, CSB-EP520564DOM-100, CSB-EP520564DOM-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	D4B2G3
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MAAQAAPELA KLDLNKNTGS VEANAVSAGG SEKEEAENEG DSDDDRDDEQ AGGSAEVNAE KKKKKRPPK KKKTAKVQSS PPRIPLTTLF PNSNFPEGEI VEYLNENSYR TTNEEKRLHD RMNNDLFTLEY RQAAEIHRQV RQYAKELIK PGATLTDIAE GIEDGVRHLT GHMGLEEGDS LVAGMGFPTG LNINHCAAHY SPNAGNKVVL QHGDVMKVDF GVHINGRIVD SA

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.