

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Caulobacter sp. Glycerol-3-phosphate dehydrogenase [NAD(P)+], Escherichia Coli Preis auf Anfrage CSB-EP531537DQR

Article Name	Recombinant Caulobacter sp. Glycerol-3-phosphate dehydrogenase [NAD(P)+], Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP531537DQR
Supplier Catalog Number	CSB-EP531537DQR
Alternative Catalog Number	CSB-EP531537DQR-1, CSB-EP531537DQR-100, CSB-EP531537DQR-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	B0T285
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MTFQHAGVIG AGAWGTALAQ VCARAGLTVT LQAREPEVTA SIAAERENTL FLAGVKLEDA ITASTDMAAL AGCDLILAVP PAQHLRSALT ALKPHLKPGA PVVLCCKGVE QGSLKLMSDV AAEVLPGAPI AVLSGPSFAG EVARGLP TAV TLACEDEALG RRIQQAALP TFRPYWADDV IGAEAGGAVK NVLAIACGIV EGKGLGRSAH AALITRGFAE LTRLAVALGG KS

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--