

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Caulobacter sp. Proline--tRNA ligase(proS),Escherichia Coli Preis auf Anfrage CSB-EP531528DQR

Article Name	Recombinant Caulobacter sp. Proline--tRNA ligase(proS),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP531528DQR
Supplier Catalog Number	CSB-EP531528DQR
Alternative Catalog Number	CSB-EP531528DQR-1, CSB-EP531528DQR-100, CSB-EP531528DQR-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	B0SZ28
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MRLSRYFLPV LKEAPSDAQI VSHQLMLRAG MIRQEAAGIY AWLPLGLRVL KKVERIVREE MDRAGAIELL MPTIQLADLW RESGRYDDYG DEMLRITDRH KRDLLYGPTA EEVVTDFIRA SVKSYKDLPK NLYNIQWKFR DERRPRFGVM RGREFLMKDA YSFDLDAEGA RKAYNRMFVA YLNLFSRMGL KAVPMRADTG PIGGDLSEHF IVLAETGESA VFCHKDLVEM AA

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--