

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Chlamydia trachomatis serovar L2 UDP-N-acetylenolpyruvoylglucosamine reductase(murB),Escherichia Coli Preis auf Anfrage CSB-EP536928DSG

Article Name	Recombinant Chlamydia trachomatis serovar L2 UDP-N-acetylenolpyruvoylglucosamine reductase(murB),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP536928DSG
Supplier Catalog Number	CSB-EP536928DSG
Alternative Catalog Number	CSB-EP536928DSG-1, CSB-EP536928DSG-100, CSB-EP536928DSG-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	B0B960
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder

Sequence	MTDSFPFSVQ ESVPLSRFST FRIGGPARYF KELTSLSEAL TVFSYLHTHP LPYIIIGKGS NCLFDDQGFD GLVLYNNIQG QEFLSDTQIK VLSGSSFALL GKRLSSQGFS GLEFAVGIPG TVGGAVFMNA GTTLANTASS LINVEIIDHS GILLSIPREK LLFSYRTSPF QKKPAFIASA TFQLTKDPQA AKRAKALIEE RILKQPYEYP SAGCIFRNPE GLSAGALIDR AG
Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.