

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Thermoanaerobacter pseudethanolicus 4-diphosphocytidyl-2-C-methyl-D-erythritol kinase(ispE),Escherichia Coli Preis auf Anfrage CSB-EP544571TMI

Article Name	Recombinant Thermoanaerobacter pseudethanolicus 4-diphosphocytidyl-2-C-methyl-D-erythritol kinase(ispE),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP544571TMI
Supplier Catalog Number	CSB-EP544571TMI
Alternative Catalog Number	CSB-EP544571TMI-1, CSB-EP544571TMI-100, CSB-EP544571TMI-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	B0KBH0
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MKKIKVKAYA KINLSLDVLG KREDGYHEIS TIMQSIDLAD ILEFEKSEIV KVFCNDHRVP LGEDNLIVKV INLLKEYQM EEGVLVKLDK RIPLAAGLAG GSADAAATIV ALDKLWNLNM SKEEKKEIAL KVGADVPFCL EGGTKLAKGI GEIFEDLKIT SMNLLLVKPD IEISTKEIYD KWDRLNFKSH HATVPVVQAI QEGNIYKIAE NIKNDLELVT SQKYGIINKI KE

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--