

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Granulibacter bethesdensis UDP-N-acetylglucosamine--N-acetylmuramyl-(pentapeptide) pyrophosphoryl-undecaprenol N-acetylglucosamine transferase(murG),Escherichia Coli Preis auf Anfrage CSB-EP603175GAAD

Article Name	Recombinant Granulibacter bethesdensis UDP-N-acetylglucosamine--N-acetylmuramyl-(pentapeptide) pyrophosphoryl-undecaprenol N-acetylglucosamine transferase(murG),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP603175GAAD
Supplier Catalog Number	CSB-EP603175GAAD
Alternative Catalog Number	CSB-EP603175GAAD-1, CSB-EP603175GAAD-100, CSB-EP603175GAAD-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	Q0BV25
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder

Sequence	MRRDQAQRPI IIAAGGTGGH FFPAEALAAE LKRRGRQIVL MTDARSGGLT STVFADTDRF VLPGAGIAGR GIRRAGQAVI ALGHGVVKAG ALLRRLEAGC IVGFGGYPCI PPVLGARLRA RNVPVILHEQ NAVLGRANRL LARKIHCLAL SFASTSHVPA GTRTLVTGNP VRPAIAALAG TSYTPPGTEG VINLLILGGS LGARVLSDVV PAALALLPPA LRQRMVRTQQ CR
Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.