

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Shewanella putrefaciens UPF0345 protein Sputcn32_0379 (Sputcn32_0379),Escherichia Coli Preis auf Anfrage CSB-EP397172STQ

Article Name	Recombinant Shewanella putrefaciens UPF0345 protein Sputcn32_0379 (Sputcn32_0379),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP397172STQ
Supplier Catalog Number	CSB-EP397172STQ
Alternative Catalog Number	CSB-EP397172STQ-1, CSB-EP397172STQ-100, CSB-EP397172STQ-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A4Y2C8
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MSLLEQVSVS KKANIYFDGK VASRSVFLAD GSKQTLGVVL PGEYEFSTSQ GEVMEVTSGR FEVLLPESTL WQEFSEGTQF ELAANVSFKI RNTAIAEYCC SYL

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.