

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Vibrio cholerae serotype O1 Vitamin B12 import ATP-binding protein BtuD(btuD),Escherichia Coli Preis auf Anfrage CSB-EP399630VEY

Article Name	Recombinant Vibrio cholerae serotype O1 Vitamin B12 import ATP-binding protein BtuD(btuD),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP399630VEY
Supplier Catalog Number	CSB-EP399630VEY
Alternative Catalog Number	CSB-EP399630VEY-1, CSB-EP399630VEY-100, CSB-EP399630VEY-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A5F1V0
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MIRVNSLQVD S RLLPLSLQA NAGEVWHVIG PNGCGKSTLL AALAGMIPFS GSVQVGGLNV SQASLSELR HRAYLAQNDK PSFQLHVFQY LALSV PANVA LERSEVASEI DQISRLLNID DKLHRSIHQL SGGEWQVRVL AGSCLQVSPV LNPSARLLIW DEPAAPLDIA QESLLYRLIE RMAGQGLTVI MANHDLNRTL RHADQVLLLS RGVLYRAGSA KEVLTQEVLQ SV

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.