

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Vibrio cholerae serotype O1 Formamidopyrimidine-DNA glycosylase(mutM),Escherichia Coli Preis auf Anfrage CSB-EP399650VEY

Article Name	Recombinant Vibrio cholerae serotype O1 Formamidopyrimidine-DNA glycosylase(mutM),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP399650VEY
Supplier Catalog Number	CSB-EP399650VEY
Alternative Catalog Number	CSB-EP399650VEY-1, CSB-EP399650VEY-100, CSB-EP399650VEY-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A5F407
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	PELPEVEVS RLGISPHLVG GTIQSLVLR PKLRWPIPQE LKQLEGQTIL AIHRRAKYLI IETAVGSAIV HLGMSGSLRI LDGDFPAKH DHVDLVM TSG KRLRYNDPRR FGAWLWCAPD ESHEVLGRLG PEPLTEAFNA EYMDKARNK RIAVKAFIMD NAAVVGVGNI YANESLFTSR LHPLRPAHSL SLEEWQTLVA NIKQVLQVAI KQGGTTLKDF TQSDGKPGYF AQE

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.