

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Recombinant Dichelobacter nodosus ATP synthase epsilon chain(atpC), Escherichia Coli Preis auf Anfrage CSB-EP399903DHZ

Article Name	Recombinant Dichelobacter nodosus ATP synthase epsilon chain(atpC),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP399903DHZ
Supplier Catalog Number	CSB-EP399903DHZ
Alternative Catalog Number	CSB-EP399903DHZ-1, CSB-EP399903DHZ-100, CSB-EP399903DHZ-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A5EXL3
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MSTPFQVYIV SPDRALFSGE VVQIVAPAEY GELGILSNHI PLIATLKPGQ VRLTKSDGEE EVLYVSGGFI EVQSKQTIIL ADEAARAAEL DEEKIKEAKA RAEKLIKSPD GEINYERMKN YEKMQIELVQ LTAQLAAIRR HRQ

## **Application Notes**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.