

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant *Flavobacterium johnsoniae* NADH-quinone oxidoreductase subunit D(nuoD), *Escherichia Coli* Preis auf Anfrage CSB-EP400017FDT

Article Name	Recombinant <i>Flavobacterium johnsoniae</i> NADH-quinone oxidoreductase subunit D(nuoD), <i>Escherichia Coli</i> Preis auf Anfrage
Biozol Catalog Number	CSB-EP400017FDT
Supplier Catalog Number	CSB-EP400017FDT
Alternative Catalog Number	CSB-EP400017FDT-1, CSB-EP400017FDT-100, CSB-EP400017FDT-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A5FKI7
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MSELLLPPEH RYAKIIKERH NEDGSELSVL NLGPTHPATH GIFQNILLMD GERILEAEPT IGYIHRAFEK IAENRPFYQI TPLTDRMNYC SSPINMGWW MTLEKLLGIE VPKRAQYLRV IVMEJARITD HIICNSILGV DTGAYTGFLY VFQFREKIYE IYEEICGARL TTNMGRIGGF ERDWSPEAFR KLDKFLEEF PAWKEFENLF ERNRIFLDRT VNVGGITAEK AM

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.