

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Recombinant Geobacter uraniireducens UPF0102 protein Gura\_3756 (Gura\_3756),Escherichia Coli Preis auf Anfrage CSB-EP400084GEI

Article Name	Recombinant Geobacter uraniireducens UPF0102 protein Gura_3756 (Gura_3756),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP400084GEI
Supplier Catalog Number	CSB-EP400084GEI
Alternative Catalog Number	CSB-EP400084GEI-1, CSB-EP400084GEI-100, CSB-EP400084GEI-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A5G7Z1
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MRSEDRADNK SLGEQGEAIA VCYLKGRKYA IVERNFRCKC GEVDIIARDG KTIVFVEVKT RRNEACGPPQ ASVTPFKQRQ ISKAALTWLA KKMLLDAPAR FDVVAILQRD HAVPEIEHIK DAFELAY

**Application Notes** 

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.