

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Geobacter uraniireducens 3-phosphoshikimate 1-carboxyvinyltransferase(aroA),Escherichia Coli Preis auf Anfrage CSB-EP402601GEI

Article Name	Recombinant Geobacter uraniireducens 3-phosphoshikimate 1-carboxyvinyltransferase(aroA),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP402601GEI
Supplier Catalog Number	CSB-EP402601GEI
Alternative Catalog Number	CSB-EP402601GEI-1, CSB-EP402601GEI-100, CSB-EP402601GEI-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A5GE08
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MQTYTIKPAQ SVRGEIAVPG DKSISHRSIM LGSIARGTTT VCGFLRGEDN MATLNAFRAM GVDVRDDGET LRIEGKGLAG LSEPTDVLDC GNSGTSMRLM TGLLAGQRFF SVLTGDQYLR NRPMKRVLEP LSRMGATIFG RAGGDKAPLA VVGRGLTGVT YASPVSSAQV KSAILLAGLY ADGETEVTEP HLSRDHSERM LRYFGADIET FAGGVKVRGR ELTGRDITVP GD

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.