

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Synechococcus sp. UPF0367 protein SynRCC307\_0258 (SynRCC307\_0258), Escherichia Coli Preis auf Anfrage CSB-EP402637SVB

|                            |                                                                                                                   |
|----------------------------|-------------------------------------------------------------------------------------------------------------------|
| Article Name               | Recombinant Synechococcus sp. UPF0367 protein SynRCC307_0258 (SynRCC307_0258), Escherichia Coli Preis auf Anfrage |
| Biozol Catalog Number      | CSB-EP402637SVB                                                                                                   |
| Supplier Catalog Number    | CSB-EP402637SVB                                                                                                   |
| Alternative Catalog Number | CSB-EP402637SVB-1, CSB-EP402637SVB-100, CSB-EP402637SVB-20                                                        |
| Manufacturer               | Cusabio                                                                                                           |
| Category                   | Proteine/Peptide                                                                                                  |
| Species Reactivity         | Bacteria                                                                                                          |
| Tag                        | The tag type will be determined during production process.                                                        |
| UniProt                    | <a href="#">A5GQK2</a>                                                                                            |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                       |
| Source                     | E.coli                                                                                                            |
| Purity                     | Greater than 85% as determined by SDS-PAGE.                                                                       |
| Form                       | Lyophilized powder                                                                                                |
| Sequence                   | MFVVELSLKL SPMPISVQRK SLEAAQGLYG EIRQAMESGH PQLMDLRCEK<br>SEEKQICLRS SEIVSVQLYE KSAMGAGSKR PGFSTGG                |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|