

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Brucella ovis Trigger factor(tig),Escherichia Coli Preis auf Anfrage CSB-EP403124BPI

Article Name	Recombinant Brucella ovis Trigger factor(tig),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP403124BPI
Supplier Catalog Number	CSB-EP403124BPI
Alternative Catalog Number	CSB-EP403124BPI-1, CSB-EP403124BPI-100, CSB-EP403124BPI-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A5VQ80
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MQVTETLNEG LKREIKVVVP AGDLEAKLAE RLETARGRAR INGFRPGKVP TAHLRKMVGK SFMAEIVNEI LNDSSRSILA ERNEKSATQP EVIMSEDEKE AEKVLGDKAD VFVSLNYEVL PAIEVKDFSK IAVTREVVDI SDEEVDEQVK RIASSTRFE TKGKAENED RVTIDYLGKL DGEPFEGGAD NDAQVLVGSG QFIPGFEEQL IGLKAGDEKV ITVTFPAEYG AA

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--