

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Shewanella sp. Queuine tRNA- ribosyltransferase(tgt),Escherichia Coli Preis auf Anfrage CSB-EP612384SAAN

|                            |                                                                                                                                                                                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Recombinant Shewanella sp. Queuine tRNA-<br>ribosyltransferase(tgt),Escherichia Coli Preis auf Anfrage                                                                                                                                                                        |
| Biozol Catalog Number      | CSB-EP612384SAAN                                                                                                                                                                                                                                                              |
| Supplier Catalog Number    | CSB-EP612384SAAN                                                                                                                                                                                                                                                              |
| Alternative Catalog Number | CSB-EP612384SAAN-1, CSB-EP612384SAAN-100, CSB-<br>EP612384SAAN-20                                                                                                                                                                                                             |
| Manufacturer               | Cusabio                                                                                                                                                                                                                                                                       |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                                              |
| Species Reactivity         | Bacteria                                                                                                                                                                                                                                                                      |
| Tag                        | The tag type will be determined during production process.                                                                                                                                                                                                                    |
| UniProt                    | <a href="#">Q0HKF7</a>                                                                                                                                                                                                                                                        |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                   |
| Source                     | E.coli                                                                                                                                                                                                                                                                        |
| Purity                     | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                   |
| Form                       | Lyophilized powder                                                                                                                                                                                                                                                            |
| Sequence                   | MKFELDTTDG RARRGRLIFE RGTVETPAFM PVGTYGTVKG MTPEEVRATG<br>ADILLGNTFH LWLRPGEEIM RKHGDLHDFM NWQRPILTDS<br>GGFQVFS LGD IRKITEEGVH FRSPINGEKI FLDPEKSMQI QHALGSDVVM<br>IFDECTPYPA TEDEARKSMQ MSLRWAKRSR DEFDRLENPN SLFGIIQGSV<br>YEDLRDESLK GLVEIGFDGY AVGGGLAVGEP KEDMHRILEH VC |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|