

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Mycobacterium tuberculosis 1D-myo-inositol 2-acetamido-2-deoxy-alpha-D-glucopyranoside deacetylase(mshB),Escherichia Coli Preis auf Anfrage CSB-EP404517MON

Article Name	Recombinant Mycobacterium tuberculosis 1D-myo-inositol 2-acetamido-2-deoxy-alpha-D-glucopyranoside deacetylase(mshB),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP404517MON
Supplier Catalog Number	CSB-EP404517MON
Alternative Catalog Number	CSB-EP404517MON-1, CSB-EP404517MON-100, CSB-EP404517MON-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A5U1L9
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder

Sequence	MSETPRLLFV HAHPDDESLS NGATIAHYTS RGAQVHVVTC TLGEEGEVIG DRWAQLTADH ADQLGGYRIG ELTAALRALG VSAPIYLGGA GRWRDSGMAG TDQRSQRRFV DADPRQTVGA LVAIIRELRP HVVVTYDPNG GYGHPDHVHT HTVTTAAVAA AGVGSGTADH PGDPWTVPKF YWTVLGLSAL ISGARALVPD DLRPEWVLPR ADEIAFGYSD DGIDAVVEAD EQ
Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.