

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Klebsiella pneumoniae subsp. pneumoniae (3R)-hydroxymyristoyl-[acyl-carrier-protein] dehydratase, Escherichia Coli Preis auf Anfrage CSB-EP408445KAX

Article Name	Recombinant Klebsiella pneumoniae subsp. pneumoniae (3R)-hydroxymyristoyl-[acyl-carrier-protein] dehydratase, Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP408445KAX
Supplier Catalog Number	CSB-EP408445KAX
Alternative Catalog Number	CSB-EP408445KAX-1, CSB-EP408445KAX-100, CSB-EP408445KAX-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6T4Y2
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MTTDTHTLHI EEILELLPHR YPFLLVDRV LDFEEGRFLRA VKNVSVNEPF FQGHFPGKPI LPGVLILEAM AQATGILAFK SVGKLEPGEL YYFAGIDEAR FKRPVVP GDQ MIM EVTFEKT RRGLTRFKGV ALVDGKVVCE ATMMCAR SRE A

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--