

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

**Recombinant Arabidopsis thaliana ATP-dependent Clp protease
proteolytic subunit-related protein 4, chloroplastic(CLPR4),Escherichia
Coli Preis auf Anfrage
CSB-EP813170DOA**

Article Name	Recombinant Arabidopsis thaliana ATP-dependent Clp protease proteolytic subunit-related protein 4, chloroplastic(CLPR4),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP813170DOA
Supplier Catalog Number	CSB-EP813170DOA
Alternative Catalog Number	CSB-EP813170DOA-1, CSB-EP813170DOA-100, CSB-EP813170DOA-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Plant
Tag	The tag type will be determined during production process.
UniProt	Q8LB10
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder

Sequence	AE SLNPLNFSSS KPKGVTMV IPFSKGS AHE QPPDLASYL FKNRIVYLGM SLVPSVTELI LAEFLYLQYE DEEKPIYLYI NSTGTTKNGE KLG YDTEAFA IYDVMGYVKP PFTLCVGNA WGEAALLLTA GAKGNRSALP SSTIMIKQPI ARFQGQATDV EIARKEIKHI KTEMVKLYSK HIGKSPEQIE ADMKRPKYFS PTEAVEYGII DKVVYNERGS QDRGVVSDLK
Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.