

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Salmonella typhimurium NH(3)-dependent NAD(+) synthetase(nadE),Escherichia Coli Preis auf Anfrage CSB-EP820798SXB

Article Name	Recombinant Salmonella typhimurium NH(3)-dependent NAD(+) synthetase(nadE),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP820798SXB
Supplier Catalog Number	CSB-EP820798SXB
Alternative Catalog Number	CSB-EP820798SXB-1, CSB-EP820798SXB-100, CSB-EP820798SXB-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	Q8ZPU5
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MTLQQEIIQA LGAKPHINPE EEIRRSVDFL KAYLKTYPF L KSLVLGISGG QDSTLAGKLS QMAIAELREE TGDNALQFIA VRLPYGVQAD EQDCQDAIAF IQPDRVLTVN IKGAVLASEQ ALREAGIELS DFVRGNEKAR ERMKAQYSIA GMTHGVVVG T DHAAEAITGF FTKYGDGGTD INPLHRLNKR QGKQLLAALG CPEHLYKKVP TADLEDDRPS LPDEAALGVT YD

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--