

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Syntrophus aciditrophicus Biotin synthase(bioB),Escherichia Coli Preis auf Anfrage CSB-EP644800SAAZ

|                            |                                                                                                                                                                                                                                                                                 |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Recombinant Syntrophus aciditrophicus Biotin synthase(bioB),Escherichia Coli Preis auf Anfrage                                                                                                                                                                                  |
| Biozol Catalog Number      | CSB-EP644800SAAZ                                                                                                                                                                                                                                                                |
| Supplier Catalog Number    | CSB-EP644800SAAZ                                                                                                                                                                                                                                                                |
| Alternative Catalog Number | CSB-EP644800SAAZ-1, CSB-EP644800SAAZ-100, CSB-EP644800SAAZ-20                                                                                                                                                                                                                   |
| Manufacturer               | Cusabio                                                                                                                                                                                                                                                                         |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                                                |
| Species Reactivity         | Bacteria                                                                                                                                                                                                                                                                        |
| Tag                        | The tag type will be determined during production process.                                                                                                                                                                                                                      |
| UniProt                    | <a href="#">Q2LY09</a>                                                                                                                                                                                                                                                          |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                     |
| Source                     | E.coli                                                                                                                                                                                                                                                                          |
| Purity                     | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                     |
| Form                       | Lyophilized powder                                                                                                                                                                                                                                                              |
| Sequence                   | MILDSTVNML AQKILDEGPQ ALSFQEACSL TELPARSIPD LFLAACKIRQ<br>HYQKDDVVL C SILNAKSGAC AENCAFCSSQ GHHR TDVVKR PLLSAEAMIE<br>EARRLDSAGA TRFSMVTSGS RLNA AEIETV GQAASAITKS TRLTVCASLG<br>QLTASSAERL KAGGVSVYHH NLETAR SFFP QICTTHDYEE DIETLLIARK<br>AGLRLCSGGI FGLGESWEQR VEMAFTLREL DV |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|