

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Solanum tuberosum Patatin group D-2, Escherichia Coli Preis auf Anfrage CSB-EP644836FIG

Article Name	Recombinant Solanum tuberosum Patatin group D-2, Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP644836FIG
Supplier Catalog Number	CSB-EP644836FIG
Alternative Catalog Number	CSB-EP644836FIG-1, CSB-EP644836FIG-100, CSB-EP644836FIG-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Plant
Tag	The tag type will be determined during production process.
UniProt	Q2MY56
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	SLEEMVT VLSIDGGGIK GIIPGTILEF LEGQLQKMDN NADARLADYF DVIGGTSTGG LLTSMITTPN ENNRPF AAAN EIVPFFFEHG PHIFNSSTGQ FFGPKYDGKY LMQVLQENLG ETRVHQALTE VAISSFDIKT NKPVIFTKSD LAKSPELDAK MYDICYSTAA APTYFP PHYF TTNTINGDKY EFNLV DGAVA TVADPALLSI SVATRLAEKD PAFASIRSLN YKKML

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--