

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Anaplasma phagocytophilum 30S ribosomal protein S20(rpsT),Escherichia Coli Preis auf Anfrage CSB-EP646198AAAQ

| Article Name | Recombinant Anaplasma phagocytophilum 30S ribosomal protein S20(rpsT),Escherichia Coli Preis auf Anfrage |
|----------------------------|--|
| Biozol Catalog Number | CSB-EP646198AAAQ |
| Supplier Catalog Number | CSB-EP646198AAAQ |
| Alternative Catalog Number | CSB-EP646198AAAQ-1, CSB-EP646198AAAQ-100, CSB-EP646198AAAQ-20 |
| Manufacturer | Cusabio |
| Category | Proteine/Peptide |
| Species Reactivity | Bacteria |
| Tag | The tag type will be determined during production process. |
| UniProt | Q2GIC8 |
| Buffer | Tris/PBS-based buffer, 6% Trehalose, pH 8.0 |
| Source | E.coli |
| Purity | Greater than 85% as determined by SDS-PAGE. |
| Form | Lyophilized powder |
| Sequence | MPNHASAKKM VRVIGKRTLS NRVRKSRVRN SVKAFVAALE ANSPIEDAIA AFRKAESNIH KCVNKGVFHR NTAARKISAL SGKLKAYDLA RLQNEQ |

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.