

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Cupriavidus pinatubonensis 23S rRNA (uracil(1939)-C(5))-methyltransferase RlmD(rlmD), Escherichia Coli Preis auf Anfrage CSB-EP670610DZR

Article Name	Recombinant Cupriavidus pinatubonensis 23S rRNA (uracil(1939)-C(5))-methyltransferase RlmD(rlmD), Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP670610DZR
Supplier Catalog Number	CSB-EP670610DZR
Alternative Catalog Number	CSB-EP670610DZR-1, CSB-EP670610DZR-100, CSB-EP670610DZR-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q46ZH7
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MPAAPRAASA SPAAVAPGSD PVVSIDSLDM EARGVGRLVN EDGTPGKVI VEGALPGETV SYKSFRKPS YEQAHLLEVR RESVMRVKPG CKHFGVCGGC SMQHLSRAQ LAIKQRVLED NLWHLISKVRP EVMFRPIAGP DWGYRYRARL TVRYVAKKGG VLVGFHERKS SYVADMTSCE ILPPHISAML VPLRELV TGL SIRD RMPQIE LAVGHEVTAL VLRILEPLTD AD

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.