

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant *Bacillus halodurans* Urease accessory protein UreD(ureD),*Escherichia Coli* Preis auf Anfrage CSB-EP873079BQS

Article Name	Recombinant <i>Bacillus halodurans</i> Urease accessory protein UreD(ureD), <i>Escherichia Coli</i> Preis auf Anfrage
Biozol Catalog Number	CSB-EP873079BQS
Supplier Catalog Number	CSB-EP873079BQS
Alternative Catalog Number	CSB-EP873079BQS-1, CSB-EP873079BQS-100, CSB-EP873079BQS-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	Q9KG55
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	<i>E.coli</i>
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MRASSHLPQL HGRLDLTFER RRGTTRLVAS EQTPPLNVSR VLRTTEEVDLA TVYLVETSGG VVSGDSQTIA IYVGEGARVE LIPQSATKVY PSRKQGESST QQVWLRVDEG AAAFWKPESI IPFRKASFLQ STRFHLHSSS TFFYGDMLTP GRVHHNERFQ YEQVDSLVEI YLDETLVVHD RVHLQPKDQL QTIGRLDGYS YYGAVWMHAP ILKGGAHLEQW MENGGQNDRF AY

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.