

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Cyanidium caldarium Imidazole glycerol phosphate synthase subunit hisH(hisH), Escherichia Coli Preis auf Anfrage CSB-EP883082DZV

Article Name	Recombinant Cyanidium caldarium Imidazole glycerol phosphate synthase subunit hisH(hisH), Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP883082DZV
Supplier Catalog Number	CSB-EP883082DZV
Alternative Catalog Number	CSB-EP883082DZV-1, CSB-EP883082DZV-100, CSB-EP883082DZV-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q9TLQ8
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MKNIGIINCG MGNLHSVSNA ISNIGFNPVI INASKDLVSF ACSALVLPGV GSFDLAIDRL EKKNLIEPVK LWIQEDRPFV GICLGLQLLF EGSDEGSKPG LKIFNGYVSQ FKHSLVKKVP HMGWKNLYFN RFNTIDNLVP HYFHYDLQPW AYFVHSYYIE PKDCYTFNTT SLTFYGKQKI VSSISYNNIL ATQFHPEKSG LFGLFILKRF LSSY

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.