

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Vibrio cholerae serotype O1
Phosphoribosylformylglycinamidine cyclo-ligase(purM),Escherichia Coli
Preis auf Anfrage
CSB-EP885196VEX

Article Name	Recombinant Vibrio cholerae serotype O1 Phosphoribosylformylglycinamidine cyclo-ligase(purM),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP885196VEX
Supplier Catalog Number	CSB-EP885196VEX
Alternative Catalog Number	CSB-EP885196VEX-1, CSB-EP885196VEX-100, CSB-EP885196VEX-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	Q9KPY6
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MSGNNPSLSY KDAGVDIDAG NALVERIKGA VKRTRRPEVM GGLGGFGALC ELPTKYKHPV LVSGTDGVGT KLRLALDMKK HDTIGIDLVA MCVNDLIVQG AEPLFFLDYY ATGKLDVDTA AEVISGIADG CLQAGCALIG GETAEMPGMY EGEDYDVAGF CVGVVEKEEI IDGSKVQVGD ALIAVGSSGP HSNGYSLVRK ILEVSKADKN ERLAGKTIGE HLLAPTKIYI KS

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--