

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Physcomitrella patens subsp. patens 50S ribosomal protein L16, chloroplastic(rpl16),Escherichia Coli Preis auf Anfrage CSB-EP751117EVD

Article Name	Recombinant Physcomitrella patens subsp. patens 50S ribosomal protein L16, chloroplastic(rpl16),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP751117EVD
Supplier Catalog Number	CSB-EP751117EVD
Alternative Catalog Number	CSB-EP751117EVD-1, CSB-EP751117EVD-100, CSB-EP751117EVD-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q6YXK9
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MLSPKRTKFR KQHRGRMKGI ATRGNSIAFG KFALQALEPS WITSRQIEAG RRAITRYARR GGKLWIRIFP DKPITMRPAE TRMGSGKGSP EYWVSVVKPG RILYEISGVP ESVARAAMRI AAYKMPIRTQ FITAITSV

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.