

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Neurospora crassa Histone acetyltransferase type B subunit 2(hat-2),Escherichia Coli Preis auf Anfrage CSB-EP752556NHA

|                            |                                                                                                                                                                                                                                                                            |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Recombinant Neurospora crassa Histone acetyltransferase type B subunit 2(hat-2),Escherichia Coli Preis auf Anfrage                                                                                                                                                         |
| Biozol Catalog Number      | CSB-EP752556NHA                                                                                                                                                                                                                                                            |
| Supplier Catalog Number    | CSB-EP752556NHA                                                                                                                                                                                                                                                            |
| Alternative Catalog Number | CSB-EP752556NHA-1, CSB-EP752556NHA-100, CSB-EP752556NHA-20                                                                                                                                                                                                                 |
| Manufacturer               | Cusabio                                                                                                                                                                                                                                                                    |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                                           |
| Species Reactivity         | Bacteria                                                                                                                                                                                                                                                                   |
| Tag                        | The tag type will be determined during production process.                                                                                                                                                                                                                 |
| UniProt                    | <a href="#">Q7S7N3</a>                                                                                                                                                                                                                                                     |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                |
| Source                     | E.coli                                                                                                                                                                                                                                                                     |
| Purity                     | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                |
| Form                       | Lyophilized powder                                                                                                                                                                                                                                                         |
| Sequence                   | MARDEIVDDV DVNMEEDDAE AEQRLINEEY KIWKKNSPFL YDMMLSTALE<br>WPTLTTQWFP DVKNPKDKSH TVHRLLLGTH TAEGKPNYLQ IAEVEIPKMV<br>ELNPRDYDEE RGEIGGYGSK ASSGEPLCIR FKITQKIDHP GEVNKARYQP<br>QNPDIATLA VDGRVLIFDR TKHSITPSGT PSPQLELIGH KEEGFGLNWN<br>PHEEGCLVTG SEDKTVLLWD LKTYEGTSKQ LK |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|