

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human UGCG protein, His tag, Unconjugated GTX00280-PRO

| Article Name | Human UGCG protein, His tag, Unconjugated |
|----------------------------|---|
| Biozol Catalog Number | GTX00280-PRO |
| Supplier Catalog Number | GTX00280-pro |
| Alternative Catalog Number | GTX00280-PRO-10 |
| Manufacturer | GeneTex |
| Category | Proteine/Peptide |
| Application | FA |
| Species Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| NCBI | 7357 |
| UniProt | Q16739 |
| Buffer | Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose, ProClin 300. |
| Expression System | E. coli |
| Form | Lyophilized powder |
| Sequence | N-terminal His-Tag, Lys39~Leu171 (NP_003349.1) |

Application Notes

UDP Glucose Ceramide Glucosyltransferase (UGCG) is an enzyme which catalyzes the first glycosylation step in glycosphingolipid biosynthesis. It belongs to the glycosyltransferase 2 family. UGCG is widely expressed and transcription is upregulated during keratinocyte differentiation. Besides, Large Multifunctional Peptidase 2 (LMP2) has been identified as an interactor of UGCG, thus a binding ELISA assay was conducted to detect the interaction of recombinant human UGCG and recombinant human LMP2. Briefly, UGCG were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl were then transferred to LMP2-coated microtiter wells and incubated for 2h at 37C. Wells were washed with PBST and incubated for 1h with anti-UGCG pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of UGCG and LMP2 was in a dose dependent manner.