

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Trident RIPA Lysis Buffer GTX400005

Article Name	Trident RIPA Lysis Buffer
Biozol Catalog Number	GTX400005
Supplier Catalog Number	GTX400005
Alternative Catalog Number	GTX400005-100
Manufacturer	GeneTex
Category	Sonstiges
Application	WB
Product Description	<p>Radio-Immunoprecipitation assay lysis buffer</p> <p>The ingredients of RIPA lysis buffer are as follows:</p> <p>25 mM Tris pH 7.6</p> <p>150 mM NaCl</p> <p>1% Triton X-100</p> <p>1% sodium deoxycholate</p> <p>1% SDS...</p>
Form	Liquid
Application Notes	<p>This product contains 1% SDS, which is higher than what is in some other commercial RIPA lysis buffers, and should facilitate a more thorough extraction of nuclear protein. The sample may be viscous following extraction due to unsheared DNA. To decrease this viscosity, we recommend increasing the volume of GTX400005 and trying either sonication or drawing the sample back-and-forth gently (to avoid bubbles) ten times through an 18G blunt needle attached to a 1ml syringe.</p>